




## Surfactant micelles in aqueous solution: critical resolution in AFM

It's a widely spread idea that performing AFM in liquids is a rather complicated research approach. Actually, many tasks related to the investigation of molecular structures and complexes in a liquid environment claim high demands on the SPM equipment. The picture below shows an AFM scan of surfactant molecules dynamically self-assembled into parallel, hemicylindrical micelles on a hydrophobic graphite surface.



**Sample:** 10 mM aqueous solution of cetyltrimethylammonium bromide surfactant in contact with HOPG.  
**AFM equipment:** [NTEGRA Prima](#), Constant force imaging mode in liquid cell using a very low set point ("soft contact mode" – force  $\sim 100$  pN).  
**AFM probe:** [CSG01](#), 0.03 N/m  
**Author:** Hannes Schniepp, The College of William and Mary, Williamsburg, VA, USA.

Observation of the dynamically self-assembled molecular structure in this case was a challenge. The imaging mode used is the "soft contact mode" – scanning with constant force using a very low set point. The vertical force of the probe needs to be kept as close as possible to the set point, in this case  $\sim 100$  pN. If the force exceeds this value, the native molecular assembly will be disturbed (the periodicity of the micelles is  $\sim 5$  nm). For lower forces, the contrast vanishes. Note that the vertical topography modulation is only 0.1-0.2 nm, but still significantly above the noise level. This example shows also that [NTEGRA nanolaboratory](#) provides enough sensitivity and stability to study dynamic molecular structures in an aqueous environment.